*highly sensitive & specific
testing in just*

30min

for Covid-19


LAMPIGEN

Fluoremetric

SARS-CoV-2

Time Line of Coronavirus

What Is The Coronavirus ?

- Diagnostic testing is the only efficient way to know the spread of the SARS-CoV-2.
 - The median incubation time is estimated to be 5.1 days with symptoms expected to be present within 12 days of infection.
- **WHO defined molecular diagnostics as the gold standard for COVID-19 testing.**

LAMPIGEN
Fluoremetric
SARS-CoV-2

Options for Response


In order to monitor and minimize the risk of COVID-19 incidence the following options for response are suggested.

- **Monitoring framework** to rapidly detect increased transmission
- **Testing strategy**
- Contact tracing
- Non-pharmaceutical interventions to prevent a large upsurge
- Risk communication

<https://www.ecdc.europa.eu/en/covid-19-pandemic>

Testing Strategies


- Comprehensive testing of all people with COVID-19-compatible symptoms.
- Testing of persons displaying symptoms compatible with COVID-19 as part of syndrome-based surveillance systems.
- Systematic random or comprehensive testing of high-risk populations or settings, irrespective of symptoms (e.g. healthcare workers).
- Testing of close contacts
- Random testing of the general population, irrespective of symptoms (point prevalence studies).
- Testing of serological markers (sero-epidemiological studies).

Travel-Related Measures to Reduce Spread

- Comprehensive testing of all people with COVID-19-compatible symptoms.
- Screening of travelers prior to or at entry at destination
 - Health questionnaires
 - Temperature screening
 - **Immunity certificate**
- **Requirement for recent negative RT-PCR test**
 - Testing before travel
 - Testing at destination
- Management of travelers with symptoms and contact tracing
- Provision of information to travelers regarding advice on travel and accessing of services

<https://www.ecdc.europa.eu/en/covid-19-pandemic>


LAMPIGEN

Fluoremetric

SARS-CoV-2


Covid-19 Testing

Covid-19 Testing

Viral gene (nucleic acid) detection is the key method for diagnosing Covid-19

Direct Method

Virus
Isolation


Viral Gene
Detection

Antigen
Detection

Indirect Method

Serology
IgM

Serology
IgG

Sensitivity


Availability

LAMPIGEN
Fluoremetric
SARS-CoV-2

Covid-19 Testing

Overview of COVID-19 Diagnostics


Covid-19 Detection Methods

Molecular Diagnostics for Infection - PCR

PCR (polymerase chain reaction) is a standard method for molecular diagnostics since 1990.

It requires specific expensive items of equipment and highly trained analysts and upwards of 4–8 h to process.

Principle of PCR Amplification Method

Covid-19 Detection Methods


Molecular Diagnostics for Infection – Isothermal Amplification


Isothermal amplification method is led to save time in diagnosis.

Five to six different types of isothermal amplification methods have been created. Among them **LAMP (loop mediated isothermal amplification)** is one of the widely experienced one.


Isothermal
Amplification
Methods


Loop-Mediated Isothermal
Amplification (LAMP)


Whole Genome
Amplification (WGA)

Strand Displacement
Amplification (SDA)

Helicase-Dependent
Amplification (HDA)

Recombinase Polymerase
Amplification (RPA)

Nucleic Acid Sequences
Based Amplification (NASBA)

Comparison of Detection Methods

LAMP has important advantages over PCR

Table 1

Comparison between PCR and LAMP based methods of viral RNA detection.
Table adapted from Nguyen *et al* (2020).

PCR	LAMP
Bulky and cumbersome	Smaller, simpler, portable.
Specialised thermal cyclers required	Only a heat block is required *
4–8 h until result	1 h until result
Requires skilled technicians	Requires no specific skill
Requires an additional reverse transcription step	Can be performed directly on RNA
Unstable reactions prone to inhibitors requiring purification steps	Stable and inhibitors tolerated, and thus purification steps not required
Detects DNA	Detects DNA and RNA
Tested on patient samples	Less tested on patient samples

LAMP is highly

**SENSITIVE,
SPECIFIC and
RAPID**

**method to detect
SARS-CoV-2**

* Heat block is not required in Fluoremetric method

ALL IN ONE

Has swabs, VTM as well as positive and negative controls in the kit

EASY TO USE

Simple to process, easy to interpret.

RAPID

Gives positive result only in 30 minutes

ACCURATE

Detects 2 different genes* It is not a subjective test (not depends on a personal evaluation)

LAMPIGEN

Fluoremetric

SARS-CoV-2

SENSITIVE

Presence of small amount of viruses are enough to detect

PORTABLE

Its device can be moved at any PoC points

FLEXIBLE

Any number of samples from 1 to 47 can be tested in one run

* N gene and ORF1ab gene

LAMPIGEN Fluorometric / SARS-CoV-2

LAMPIGEN Fluorometric

- LAMPIGEN is a vitro **Point of Care** diagnostic test for **SARS-CoV-2** virus based on Reverse Transcription Loop-mediated isothermal Amplification (**RT-LAMP**) technology with the fluorometric detection method.
- LAMPIGEN is associated with fluorometric RT-LAMP devices (Genie® II, III and HT)

LAMPIGEN Components


LAMPIGEN Components

LAMPIGEN supplies all necessary components in its kit box

Including;

- Oro/nasopharyngeal swabs &
- Virus transport medium (VTM) units.

LAMPIGEN Components

LAMPIGEN Fluorometric SARS-CoV2 RT-LAMP Kit Components *

Components
Stored at
-15°C to -25°C

Components	Number of vials	Reactions per vial
RT-LAMP Mastermix	10	100
Primer Mix-1 (N gene)	2	250
Primer Mix-2 (ORF1ab gene)	2	250
Positive Control	1	20
Negative Control	1	20

Components
Stored at
Room temperature

Components	Number of sets	Numbers per set	Total units
Tube-Strip	128	8	1024
Virus Transport Medium (VTM)			500
Oro/Nasopharyngeal swabs			500

* 500 tests/kit


Test Process – Step 1

Sample Collection

Swab should be inserted first into one nostril of the patient up to 2.5 cm (1 inch) from the edge of the nostril.

It must be rolled 5 times along the mucosa inside the nostril to ensure that both mucus and cells are collected.

Then using the same swab process should be repeated for the other nostril to ensure that an adequate sample is collected from both nasal cavities.


Test Process – Step 2

VTM

Then swab is placed in VTM.

The tail of the swab is broken and the body is left into the VTM.

Both **swabs** and **VTM** are present in the **KIT Box**

Test Process – Step 3

Preparation of Testing Samples


Since two genes will be tested samples should be duplicated into two tubes.

First; put 17,5 µl of **Mastermix** into two adjacent tubes.

Then add 2,5 µl of **Primer-1** into the first tube and 2,5 µl of **Primer-2** into the second tube.

A fresh reaction mix should be prepared before each batch of samples is tested.

Samples should be tested in duplicate: wells 1&2; wells 3&4; wells 5&6; wells 7&8. Duplicates must be set up in the configuration shown in the figure. One pair of controls should be placed in each run.

Sensitive and Accurate Detection

PERFORMANCE

**FAST AND ACCURATE TESTING FOR
COVID -19**

Detection Sensitivity	3 - 5 copies / μ l
Genes Detected	N gene & ORF1ab gene
Detection Method	Fluorometric

LAMPIGEN detects **2 genes**
with **Fluorometric Method**
which is sensitive enough for only **3 to 5 copies in μ l**

LAMPIGEN
Fluoremetric
SARS-CoV-2

Sensitive and Accurate Detection

Highly Specific to SARS-CoV2

Tested Respiratory Pathogens		
Legionella pneumophila	Human Coronavirus HKU1	Human Parainfluenza 3
Haemophilus influenzae	Influenza A virus (Flu A)	Human Parainfluenza 4
Streptococcus pneumoniae	Influenza B virus (Flu B)	Human bocavirus
Bordetella parapertussis	Influenza A/H1	Respiratory syncytial virus
Human Coronavirus NL63	Influenza A/H3	Human adenovirus
Human Coronavirus 229E	Human Parainfluenza 1	Enterovirus
Human Coronavirus OC43	Human Parainfluenza 2	Parechovirus

All the above pathogens were tested with **LAMPIGEN** to check for cross reactions and resulted as **NEGATIVE**.

No cross reaction was observed with any of the possible pathogenic microorganisms during development tests.

High Sensitivity & High Specificity

LAMPIGEN has high sensitivity & high specificity

Overall Diagnostic Sensitivity ¹	100 %
Overall Diagnostic Specificity ²	97 %

1- 100% (samples with <CT25) (Concordance with RT-qPCR to detect positives)

2- Samples were tested in singles using the qRT-PCR. For the statistics, it was assumed that the qRT-PCR results were correct. Cycle threshold (CT) is the number of cycles required for the fluorescent signal to cross a defined threshold for the comparator qRT-PCR.

(Concordance with RT-qPCR to detect negatives)

LAMPIGEN RT LAMP PCRCOVID-19 Kit

Qualitative detection of SARS-CoV-2 RNA

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CE

LAMPIGEN
SARS-CoV-2

ATTESTATION OF COMPLIANCE UYGUNLUK ONAYI

Üretici / Manufacturer
PHARMALINE SAĞLIK HİZMETLERİ TİCARET A.Ş.
Adres / Address
SAHRAYICEDİT MAH. ATATÜRK CAD. NO: 43/6 KADIKÖY / İSTANBUL / TÜRKİYE
Telefon / Phone
+90 216 345 86 66
E-mail / Web
info@pharmaline.com.tr / www.pharmaline.com.tr

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The following described product in our delivered version complies with appropriate basic safety and health requirements of 98/79/EEC In Vitro Diagnostic Medical Devices Directive on its design and type, as brought into circulation by us, in case of alteration of the product, not agreed upon by us, this declaration will lose its validity.

Ürün İsmi- Product Name
LAMPİGEN COVID-19 RT-LAMP PCR KIT
LAMPİGEN VIRAL TRANSPORT MEDIUM KIT

COVIDİGEN COVID-19 QRT PCR KIT
COVIDİGEN VIRAL TRANSPORT MEDIUM KIT

Modeller / Models
LAMPİGEN LAMPPH2020-500, PH2021VTM
COVIDİGEN PH2020-100, PH2020-500, PH2020VTM
Ürünün Markası / Product Brand
LAMPİGEN / COVIDİGEN

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Accredited System Certification Approval

Seal and Date
28.07.2020

Certification Manager
MRUDULA BABAU

AB-UYGUNLUK BEYANI EC-DECLARATION OF CONFORMITY

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ÜRETİCİNİN ADRESİ / ADDRESS OF MANUFACTURER:

ATATÜRK CADDESİ, ATATÜRK CADDESİ KÖNÜTLARİ NO:43 A BLOK DAİRE:6
SAHRAYICEDİT MAHALLESİ -KADIKÖY -İSTANBUL-TÜRKİYE

ÜRÜN TANIMI / PRODUCT DESCRIPTION:

Lampigen COVID-19 RT-LAMP PCR Kit [Catalog No: LAMPPH2020-500]

Lampigen COVID-19 RT-LAMP PCR Kit [Catalog No: LAMPPH2020-500]

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98/79/EC IN VITRO DIAGNOSTIC MEDICAL DEVICES DIRECTIVE ANNEX III PRODUCT

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WE HEREBY DECLARE THAT THE DISTRIBUTED CE MARKED MEDICAL DEVICES LISTED BELOW CONFORM TO THE REQUIRED TECHNICAL DOCUMENTATION, IN ACCORDANCE WITH ANNEX III OF THE COUNCIL DIRECTIVE 98/79/EC, CONCERNING IN VITRO DIAGNOSTIC MEDICAL DEVICES.

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CLASSIFIED AS GENERAL/OTHER IVD DEVICES
THE JUSTIFICATION FOR COMPLIANCE WITH THIS DIRECTIVE IS GIVEN ON THE PRODUCT TECHNICAL FILE.

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